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Using 2D experiments with reverse heteronuclear correlation on a benchtop NMR spectrometer X-Pulse to study the structure of gemfibrozil

The article is presented by the company LLC "AVRORA"

The article discusses application of two powerful NMR techniques ME-HSQC and HMBC to determine the structure of gemfibrozil molecules as an example on the benchtop NMR spectrometer X-Pulse. The obtained results give grounds to assert that the use of the benchtop X-Pulse NMR spectrometer allows obtaining reliable results in a short time.

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¹ Pagès G., Gerdova A., Williamson D., Gilard V., Martino R., Malet-Martino M. Evaluation of a Benchtop Cryogen-Free Low-Field 1H NMR Spectrometer for the Analysis of Sexual Enhancement and Weight Loss Dietary Supplements Adulterated with Pharmaceutical Substances. *ACS Analytical Chemistry*. 2014;86(23):11897–11904.

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⁴ Caprio V., McLachlan A. S., Sutcliffe O. B., Williamson D. C., Mewis R. E. Structural Elucidation of Unknowns: A Spectroscopic Investigation with an Emphasis on 1D and 2D 1H Nuclear Magnetic Resonance Spectroscopy. *Chemistry: Bulgarian Journal of Science Education*. 2018;27(2):223–234.

⁵ Fallaise D., Tweedie H. B., Konzuk J., Cheyne C., Mack E. E., Longstaffe J. G. Practical application of 1 H benchtop NMR spectroscopy for the characterization of a nonaqueous phase liquid from a contaminated environment. *Magnetic Resonance in Chemistry*. 2019;57(2–3):93–100.

⁶ Antonides L. H., Brignall R. M., Costello A., Ellison J., Firth S. E., Gilbert N., Groom B. J., Hudson S. J., Hulme M. C., Marron J., Pullen Z. A., Robertson T. B. R., Schofield C. J., Williamson D. C., Kemsley E. K., Sutcliffe O. B., Mewis R. E. Rapid Identification of Novel Psychoactive and Other Controlled Substances Using Low-Field 1H NMR Spectroscopy. *ACS Omega*. 2019;4(4):7103–7112.

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¹ Claridge T.D.W. High-Resolution NMR Techniques in Organic Chemistry. Volume 2. 2nd Edition. Amsterdam: Elsevier Science. 2008. 398 с.

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Figure 1. HMBC spectrum for 1 M gemfibrozil solution. 16 scans with a total experiment time of 1 hour and 40 minutes. The spectrum was obtained on a desktop NMR spectrometer X-Pulse

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Figure 2. ME-HSQC spectrum for 1 M gemfibrozil solution. 4 scans with a total experiment time of 35 min. The spectrum was obtained on the desktop NMR spectrometer X-Pulse. The numbers on the cross peaks correspond to the numbered positions of the carbon atoms in figure 1

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Figure 3. HMBC spectrum for 1 M gemfibrozil solution. 16 scans with a total experiment time of 1 hour and 40 minutes. The spectrum was obtained on a desktop NMR spectrometer X-Pulse

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